

For Research Use Only

PRDX3 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 66810-2-PBS



Basic Information

Catalog Number:

66810-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1062

GenBank Accession Number:

BC007062

GeneID (NCBI):

10935

UNIPROT ID:

P30048

Full Name:

peroxiredoxin 3

Calculated MW:

27 kDa

Purification Method:

Protein G purification

CloneNo.:

1G9A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50218-1, PRDX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-2-PBS. Detection antibody: 66810-3-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50218-2, PRDX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-1-PBS. Detection antibody: 66810-2-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL.