

For Research Use Only

MYH14 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66825-4-PBS



Basic Information

Catalog Number:

66825-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28063

GenBank Accession Number:

NM_001077186

GeneID (NCBI):

79784

UNIPROT ID:

Q7Z406

Full Name:

myosin, heavy chain 14

Calculated MW:

228 kDa

Purification Method:

Protein G purification

CloneNo.:

2A9A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50505-2, MYH14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66825-4-PBS. Detection antibody: 66825-5-PBS. Standard:Ag28063. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50505-4, MYH14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66825-4-PBS. Detection antibody: 66825-7-PBS. Standard:Ag28063. Range: 0.391-100 ng/mL