

For Research Use Only

IRF5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66835-3-PBS



Basic Information

Catalog Number:

66835-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5442

GenBank Accession Number:

BC004201

GeneID (NCBI):

3663

UNIPROT ID:

Q13568

Full Name:

interferon regulatory factor 5

Calculated MW:

56 kDa

Purification Method:

Protein G purification

CloneNo.:

3D1E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51481-1, IRF5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66835-2-PBS. Detection antibody: 66835-3-PBS.



Cytometric bead array standard curve of MP51481-1, IRF5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66835-2-PBS. Detection antibody: 66835-3-PBS. Standard:Ag5442. Range: 1.563-200 ng/mL.



Bi-layer interferometry (BLI) kinetic assays of 66835-3-Ig against Human IRF5 were performed. The affinity constant is 0.915 nM.