

For Research Use Only

# NeuN Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66836-2-PBS



## Basic Information

Catalog Number:	66836-2-PBS	GenBank Accession Number:	NM_001082575	Purification Method:	Protein G purification
Source:	Mouse	GeneID (NCBI):	146713	CloneNo.:	3H2B2
Isotype:	IgG1	UNIPROT ID:	A6NFN3		
Immunogen Catalog Number:	AG28016	Full Name:	hexaribonucleotide binding protein 3		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

