

For Research Use Only

RICTOR Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66867-3-PBS



Basic Information

Catalog Number:

66867-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28487

GenBank Accession Number:

BC029608

GeneID (NCBI):

253260

UNIPROT ID:

Q6R327

Full Name:

rapamycin-insensitive companion of mTOR

Calculated MW:

192 kDa

Purification Method:

Protein G purification

CloneNo.:

1E3E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50234-1, RICTOR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66867-2-PBS. Detection antibody: 66867-3-PBS. Standard: Ag28487. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 66867-3-PBS against Human RICTOR were performed. The affinity constant is 14.1 nM.