

For Research Use Only

CD8A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66868-2-PBS

Basic Information

Catalog Number:

66868-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11111

GenBank Accession Number:

BC025715

GeneID (NCBI):

925

ENSEMBL Gene ID:

ENSG00000153563

UNIPROT ID:

P01732

Full Name:

CD8a molecule

Calculated MW:

235 aa, 26 kDa

Purification Method:

Protein G purification

CloneNo.:

1E12B3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50246-1, CD8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66868-2-PBS. Detection antibody: 66868-3-PBS. Standard:Ag11111. Range: 3.125-100 ng/mL.