

For Research Use Only

ANXA10 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66869-2-PBS



Basic Information

Catalog Number: 66869-2-PBS	GenBank Accession Number: BC007320	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 11199	CloneNo.: 2G5A2
Source: Mouse	UNIPROT ID: Q9UJ72	
Isotype: IgG1	Full Name: annexin A10	
Immunogen Catalog Number: AG27098	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51324-1, ANXA10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66869-2-PBS. Detection antibody: 66869-3-PBS.



Cytometric bead array standard curve of MP51324-1, ANXA10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66869-2-PBS. Detection antibody: 66869-3-PBS. Standard:Ag27098. Range: 0.195-200 ng/mL.