

For Research Use Only

ANXA10 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66869-4-PBS



Basic Information

Catalog Number:	66869-4-PBS	GenBank Accession Number:	BC007320	Purification Method:	Protein G Magarose purification
Concentration:	1 mg/ml	GeneID (NCBI):	11199	CloneNo.:	3B2B10
Source:	Mouse	UNIPROT ID:	Q9UJ72		
Isotype:	IgG1	Full Name:	annexin A10		
Immunogen Catalog Number:	AG27098	Calculated MW:	37 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51324-2, ANXA10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66869-4-PBS. Detection antibody: 66869-5-PBS.



Cytometric bead array standard curve of MP51324-2, ANXA10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66869-4-PBS. Detection antibody: 66869-5-PBS. Standard:Ag27098. Range: 1.563-200 ng/mL.