

For Research Use Only

# SREBF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66875-2-PBS



## Basic Information

Catalog Number:

66875-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5484

GenBank Accession Number:

BC063281

GeneID (NCBI):

6720

UNIPROT ID:

P36956

Full Name:

sterol regulatory element binding  
transcription factor 1

Calculated MW:

1177 aa, 125 kDa

Purification Method:

Protein G purification

CloneNo.:

1A9G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50853-1, SREBF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66875-2-PBS. Detection antibody: 66875-3-PBS. Standard:Ag5484. Range: 0.195-12.5 ng/mL