

For Research Use Only

ASF1A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66884-2-PBS



Basic Information

Catalog Number:

66884-2-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG1233

GenBank Accession Number:

BC010878

GeneID (NCBI):

25842

UNIPROT ID:

Q9Y294

Full Name:

ASF1 anti-silencing function 1 homolog A (S. cerevisiae)

Calculated MW:

204 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

4G1C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51470-1, ASF1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66884-2-PBS. Detection antibody: 66884-1-PBS.



Cytometric bead array standard curve of MP51470-1, ASF1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66884-2-PBS. Detection antibody: 66884-1-PBS. Standard: Ag1233. Range: 0.781-100 ng/mL.