

For Research Use Only

MTOR Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66888-2-PBS

Basic Information

Catalog Number:

66888-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28395

GenBank Accession Number:

NM_004958

GeneID (NCBI):

2475

UNIPROT ID:

P42345

Full Name:

FK506 binding protein 12-rapamycin associated protein 1

Calculated MW:

289 kDa

Purification Method:

Protein A purification

CloneNo.:

1B12D1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50209-1, MTOR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66888-2-PBS. Detection antibody: HRP-conjugated 66888-3-PBS. Standard:Ag28395. Range: 0.049-50 ng/mL.



Sandwich ELISA standard curve of MP50209-2, MTOR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66888-2-PBS. Detection antibody: HRP-conjugated 66888-4-PBS. Standard:Ag28395. Range: 0.049-50 ng/mL.