

For Research Use Only

DMXL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66891-3-PBS



Basic Information

Catalog Number:

66891-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG26302

GenBank Accession Number:

BC140781

GeneID (NCBI):

23312

UNIPROT ID:

Q8TDJ6

Full Name:

Dmx-like 2

Calculated MW:

3036 aa, 340 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F5B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP1508-1, DMXL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66891-1-PBS. Detection antibody: 66891-3-PBS.



Cytometric bead array standard curve of MP1508-1, DMXL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66891-1-PBS. Detection antibody: 66891-3-PBS. Standard:Ag26302. Range: 1.563-200 ng/mL.