

For Research Use Only

Glucocorticoid receptor Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66904-2-PBS



Basic Information

Catalog Number:	66904-2-PBS	GenBank Accession Number:	BC015610	Purification Method:	Protein G purification
Source:	Mouse	GeneID (NCBI):	2908	CloneNo.:	2E8B5
Isotype:	IgG1	UNIPROT ID:	P04150		
Immunogen Catalog Number:	AG28431	Full Name:	nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)		
		Calculated MW:	86 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS.



Cytometric bead array standard curve of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS. Standard:Ag28431. Range: 0.781-100 ng/mL.