

For Research Use Only

Glucocorticoid receptor Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66904-3-PBS



Basic Information

Catalog Number:

66904-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28431

GenBank Accession Number:

BC015610

GeneID (NCBI):

2908

UNIPROT ID:

P04150

Full Name:

nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)

Calculated MW:

86 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B10B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS.



Cytometric bead array standard curve of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS. Standard:Ag28431. Range: 0.781-100 ng/mL.