

For Research Use Only

CDC25C Monoclonal antibody, PBS Only (Detector)

Catalog Number:66912-3-PBS



Basic Information

Catalog Number:

66912-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10240

GenBank Accession Number:

BC019089

GeneID (NCBI):

995

UNIPROT ID:

P30307

Full Name:

cell division cycle 25 homolog C (S. pombe)

Calculated MW:

473 aa, 53 kDa

Purification Method:

Protein G purification

CloneNo.:

1F11C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

