

For Research Use Only

NFATC3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66924-2-PBS

Basic Information

Catalog Number: 66924-2-PBS	GenBank Accession Number: BC001050	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 4775	CloneNo.: 1C3B10
Source: Mouse	UNIPROT ID: Q12968	Affinity: $K_D = 5.73 \times 10^{-11}$
Isotype: IgG1	Full Name: nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3	$K_{off} = 4.88 \times 10^{-6}$ $K_{on} = 8.53 \times 10^4$
Immunogen Catalog Number: AG27306	Calculated MW: 116 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP50022-1, NFATC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66924-2-PBS. Detection antibody: HRP-conjugated 66924-1-PBS. Standard: Ag27306. Range: 0.625-20 ng/mL.