

For Research Use Only

Kir4.1 Monoclonal antibody, PBS Only



Catalog Number: 66931-1-PBS

Basic Information

Catalog Number:

66931-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16785

GenBank Accession Number:

BC034036

GeneID (NCBI):

3766

UNIPROT ID:

P78508

Full Name:

potassium inwardly-rectifying channel, subfamily J, member 10

Calculated MW:

379 aa, 41 kDa

Observed MW:

40 kDa

Purification Method:

Protein G purification

CloneNo.:

2D6A6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 66931-1-Ig (Kir4.1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66931-1-PBS in a different storage buffer formulation.