

For Research Use Only

LDLRAP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66932-3-PBS



Basic Information

Catalog Number:

66932-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28608

GenBank Accession Number:

BC029770

GeneID (NCBI):

26119

UNIPROT ID:

Q5SW96

Full Name:

low density lipoprotein receptor adaptor protein 1

Calculated MW:

308 aa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B2A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51596-1, LDLRAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66932-3-PBS. Detection antibody: 66932-2-PBS.



Cytometric bead array standard curve of MP51596-1, LDLRAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66932-3-PBS. Detection antibody: 66932-2-PBS. Standard:Ag28608. Range: 1.563-200 ng/mL