

For Research Use Only

RACK1; GNB2L1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66940-3-PBS



Basic Information

Catalog Number:

66940-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG26708

GenBank Accession Number:

BC019093

GeneID (NCBI):

10399

UNIPROT ID:

P63244

Full Name:

guanine nucleotide binding protein (G protein), beta polypeptide 2-like 1

Calculated MW:

36 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3D8E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

