

For Research Use Only

# DBR1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66951-2-PBS



## Basic Information

|                                     |                                                            |                             |
|-------------------------------------|------------------------------------------------------------|-----------------------------|
| Catalog Number:<br>66951-2-PBS      | GenBank Accession Number:<br>BC009472                      | Purification Method:<br>N/A |
| Source:<br>Mouse                    | GeneID (NCBI):<br>51163                                    | CloneNo.:<br>3B10B2         |
| Isotype:<br>IgG1                    | UNIPROT ID:<br>Q9UK59                                      |                             |
| Immunogen Catalog Number:<br>AG8830 | Full Name:<br>debranching enzyme homolog 1 (S. cerevisiae) |                             |
|                                     | Calculated MW:<br>544 aa, 62 kDa                           |                             |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array sample test of MP51446-1, DBR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66951-2-PBS. Detection antibody: 66951-3-PBS.



Cytometric bead array standard curve of MP51446-1, DBR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66951-2-PBS. Detection antibody: 66951-3-PBS. Standard:Ag8830. Range: 0.391-100 ng/mL.