

For Research Use Only

LXN Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66957-2-PBS

Basic Information

Catalog Number:

66957-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27922

GenBank Accession Number:

BC005346

GeneID (NCBI):

56925

UNIPROT ID:

Q9BS40

Full Name:

latexin

Calculated MW:

222 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

1G2A4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP50023-1, LXN
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 66957-1-PBS. Detection
antibody: HRP-conjugated 66957-2-PBS. Standard:
Ag27922. Range: 0.313-20 ng/mL.