

For Research Use Only

RRAS Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66959-3-PBS



Basic Information

Catalog Number: 66959-3-PBS	GenBank Accession Number: BC016318	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 6237	CloneNo.: 2A7C12
Isotype: IgG2b	UNIPROT ID: P10301	
Immunogen Catalog Number: AG26564	Full Name: related RAS viral (r-ras) oncogene homolog	
	Calculated MW: 218 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51610-1, RRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66959-2-PBS. Detection antibody: 66959-3-PBS.



Cytometric bead array standard curve of MP51610-1, RRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66959-2-PBS. Detection antibody: 66959-3-PBS. Standard:Ag26564. Range: 0.781-100 ng/mL.