

For Research Use Only

CD40 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66965-2-PBS



Basic Information

Catalog Number:

66965-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28040

GenBank Accession Number:

BC012419

GeneID (NCBI):

958

ENSEMBL Gene ID:

ENSG00000101017

UNIPROT ID:

P25942

Full Name:

CD40 molecule, TNF receptor superfamily member 5

Calculated MW:

277 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

5G1A8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50064-1, CD40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66965-2-PBS. Detection antibody: 65103-1-PBS. Standard: Ag28040. Range: 78.1-5000 pg/mL.