

For Research Use Only

IKZF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66966-2-PBS



Basic Information

Catalog Number: 66966-2-PBS	GenBank Accession Number: BC018349	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 10320	CloneNo.: 1C1G12
Isotype: IgG1	UNIPROT ID: Q13422	
Immunogen Catalog Number: AG28637	Full Name: IKAROS family zinc finger 1 (Ikaros)	
	Calculated MW: 477 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS.



Cytometric bead array standard curve of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS. Standard: Ag28637. Range: 1.563-200 ng/mL.