

For Research Use Only

# IKZF1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66966-3-PBS



## Basic Information

<b>Catalog Number:</b> 66966-3-PBS	<b>GenBank Accession Number:</b> BC018349	<b>Purification Method:</b> Protein G purification
<b>Source:</b> Mouse	<b>GeneID (NCBI):</b> 10320	<b>CloneNo.:</b> 2C9G5
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> Q13422	
<b>Immunogen Catalog Number:</b> AG28637	<b>Full Name:</b> IKAROS family zinc finger 1 (Ikaros)	
	<b>Calculated MW:</b> 477 aa, 53 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS.



Cytometric bead array standard curve of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS. Standard: Ag28637. Range: 1.563-200 ng/mL.