

Hexokinase 2 Monoclonal antibody

Catalog Number: 66974-1-Ig 56 Publications

Basic Information

Catalog Number: 66974-1-Ig	GenBank Accession Number: BC021116	Purification Method: Protein G purification
Concentration: 1500 ug/ml	GeneID (NCBI): 3099	CloneNo.: 2A11C3
Source: Mouse	UNIPROT ID: P52789	Recommended Dilutions: WB 1:5000-1:20000 IHC 1:250-1:1000 IF/ICC 1:375-1:1500
Isotype: IgG1	Full Name: hexokinase 2	
Immunogen Catalog Number: AG16895	Calculated MW: 102 kDa Observed MW: 102 kDa	

Applications

Tested Applications: WB, IHC, IF/ICC, FC (Intra), ELISA	Positive Controls: WB : HEK-293 cells, NIH/3T3 cells, LNCaP cells, HeLa cells, HepG2 cells, Jurkat cells, K-562 cells, 4T1 cells IHC : mouse skeletal muscle tissue, IF/ICC : HeLa cells,
Cited Applications: WB, IHC, IF	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse, rat, pig	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Hexokinase 2 (HK2), a rate-limiting enzyme in the first step of glycolysis pathway, expresses at high level in cancer cells compared with normal cells. HK2 provides a new target for cancer therapy due to its pivotal role in tumor tumorigenic and metastatic process. HK1 is constitutively expressed in most mammalian adult tissues. HK2, however, although is abundantly expressed in embryonic tissues, is expressed at high levels only in limited number of adult tissues such as adipose, skeletal, and cardiac muscles. (PMID: 29305912, PMID: 28427443, PMID: 23911236). The antibody also recognizes Hexokinase 1.

Notable Publications

Author	Pubmed ID	Journal	Application
Mengyuan Li	34558014	Virol Sin	WB
Yanling Gou	36069580	Immunology	IHC
Yubo Zhao	36064939	Cell Death Dis	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 66974-1-Ig (Hexokinase 2 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Lamin B1 Monoclonal antibody (HRP-66095) as loading control.



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 66974-1-Ig (Hexokinase 2 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Hexokinase 2 antibody (66974-1-Ig, Clone: 2A11C3) at dilution of 1:750 and CoraLite®488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1).



1x10⁶ HepG2 cells were intracellularly stained with 0.25 ug Hexokinase 2 Monoclonal antibody (66974-1-Ig, Clone:2A11C3) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.25 ug Mouse IgG1 Isotype Control (MOPC-21) (65124-1-Ig, Clone: MOPC-21) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).