

For Research Use Only

fetuin-B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67002-3-PBS



Basic Information

Catalog Number: 67002-3-PBS	GenBank Accession Number: BC074734	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 26998	CloneNo.: 1F3A7
Isotype: IgG1	UNIPROT ID: Q9UGM5	
Immunogen Catalog Number: AG12463	Full Name: fetuin B	
	Calculated MW: 382 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array sample test of MP51617-1, Fetuin-B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67002-2-PBS. Detection antibody: 67002-3-PBS.



Cytometric bead array sample test of MP51617-2, Fetuin-B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67002-4-PBS. Detection antibody: 67002-3-PBS.



Cytometric bead array standard curve of MP51617-1, Fetuin-B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67002-2-PBS. Detection antibody: 67002-3-PBS. Standard: Ag12463. Range: 3.125-400 ng/mL.



Cytometric bead array standard curve of MP51617-2, Fetuin-B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67002-4-PBS. Detection antibody: 67002-3-PBS. Standard: Ag12463. Range: 3.125-400 ng/mL.