

For Research Use Only

RRM2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67006-3-PBS



Basic Information

Catalog Number:

67006-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28664

GenBank Accession Number:

BC030154

GeneID (NCBI):

6241

UNIPROT ID:

P31350

Full Name:

ribonucleotide reductase M2

polypeptide

Calculated MW:

389 aa, 45 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C6B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51399-1, RRM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67006-2-PBS. Detection antibody: 67006-3-PBS.



Cytometric bead array standard curve of MP51399-1, RRM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67006-2-PBS. Detection antibody: 67006-3-PBS. Standard:Ag28664. Range: 1.563-200 ng/mL.