

For Research Use Only

# GNAI2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67007-2-PBS

## Basic Information

Catalog Number:

67007-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28560

GenBank Accession Number:

BC012138

GeneID (NCBI):

2771

UNIPROT ID:

P04899

Full Name:

guanine nucleotide binding protein (G protein), alpha inhibiting activity polypeptide 2

Calculated MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1C11D2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve ofMP50015-1, GNAI2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67007-2-PBS. Detection antibody: HRP-conjugated 67007-3-PBS. Standard: Ag28560. Range: 0.313-20 ng/mL.