

For Research Use Only

# EDN1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67008-5-PBS



## Basic Information

Catalog Number:

67008-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28334

GenBank Accession Number:

BC009720

GeneID (NCBI):

1906

UNIPROT ID:

P05305

Full Name:

endothelin 1

Calculated MW:

212 aa, 24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E4F11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sample test of MP51065-2, EDN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67008-4-PBS. Detection antibody: 67008-5-PBS.



Cytometric bead array standard curve of MP51065-2, EDN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67008-4-PBS. Detection antibody: 67008-5-PBS. Standard: Ag28334. Range: 1.563-200 ng/mL.