

For Research Use Only

COMMD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67016-2-PBS



Basic Information

Catalog Number:

67016-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28596

GenBank Accession Number:

BC022046

GeneID (NCBI):

150684

UNIPROT ID:

Q8N668

Full Name:

copper metabolism (Murr1) domain containing 1

Calculated MW:

190 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G7D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51604-1, COMMD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67016-1-PBS. Detection antibody: 67016-2-PBS.



Cytometric bead array standard curve of MP51604-1, COMMD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67016-1-PBS. Detection antibody: 67016-2-PBS. Standard:Ag28596. Range: 0.313-40 ng/mL