

For Research Use Only

# STBD1 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 67018-2-PBS

## Basic Information

Catalog Number: 67018-2-PBS	GenBank Accession Number: BC022301	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 8987	CloneNo.: 1B3B3
Source: Mouse	UNIPROT ID: O95210	
Isotype: IgG1	Full Name: starch binding domain 1	
Immunogen Catalog Number: AG28613	Calculated MW: 41 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50130-1, STBD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67018-2-PBS. Detection antibody: 67018-3-PBS. Standard:Ag28613. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50130-2, STBD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67018-1-PBS. Detection antibody: 67018-2-PBS. Standard:Ag28613. Range: 1.563-100 ng/mL