

For Research Use Only

# TIMELESS Monoclonal antibody, PBS Only (Capture)

Catalog Number:67022-2-PBS



## Basic Information

Catalog Number:

67022-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6283

GenBank Accession Number:

BC050557

GeneID (NCBI):

8914

UNIPROT ID:

Q9UNS1

Full Name:

timeless homolog (Drosophila)

Calculated MW:

139 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E5A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51685-1, TIMELESS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67022-2-PBS. Detection antibody: 67022-1-PBS.



Cytometric bead array standard curve of MP51685-1, TIMELESS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67022-2-PBS. Detection antibody: 67022-1-PBS. Standard: Ag6283. Range: 0.313-40 ng/mL.