

For Research Use Only

CUL7 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67034-2-PBS



Basic Information

Catalog Number:

67034-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6943

GenBank Accession Number:

BC033647

GeneID (NCBI):

9820

UNIPROT ID:

Q14999

Full Name:

cullin 7

Calculated MW:

1698 aa, 191 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F6B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51627-1, CUL7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67034-2-PBS. Detection antibody: 67034-3-PBS.



Cytometric bead array standard curve of MP51627-1, CUL7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67034-2-PBS. Detection antibody: 67034-3-PBS. Standard:Ag6943. Range: 0.781-100 ng/mL