

For Research Use Only

# MED30 Monoclonal antibody, PBS Only



Catalog Number: 67038-1-PBS

Featured Product

## Basic Information

Catalog Number:

67038-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG10235

GenBank Accession Number:

BC008226

GeneID (NCBI):

90390

UNIPROT ID:

Q96HR3

Full Name:

mediator complex subunit 30

Calculated MW:

178 aa, 20 kDa

Observed MW:

20-25 kDa

Purification Method:

Protein A purification

CloneNo.:

1D7C4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, rat, pig

## Background Information

MED30, also named as THRAP6 and TRAP25, is a component of the Mediator complex. MED30 has two isoforms with MW 16-20 kDa and 20-25 kDa.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67038-1-Ig (MED30 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67038-1-PBS in a different storage buffer formulation.



WB result of MED30 antibody (67038-1-Ig; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MED30 transfected HEK-293 cells. This data was developed using the same antibody clone with 67038-1-PBS in a different storage buffer formulation.