

For Research Use Only

# PELP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67050-2-PBS



## Basic Information

Catalog Number:

67050-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25729

GenBank Accession Number:

BC069058

GeneID (NCBI):

27043

UNIPROT ID:

Q8IZL8

Full Name:

proline, glutamate and leucine rich protein 1

Calculated MW:

1130 aa, 120 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D11G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data

