

For Research Use Only

PELP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67050-3-PBS



Basic Information

Catalog Number: 67050-3-PBS	GenBank Accession Number: BC069058	Purification Method: Protein A Magarose purification
Source: Mouse	GeneID (NCBI): 27043	CloneNo.: 1G8E9
Isotype: IgG2a	UNIPROT ID: Q8IZL8	
Immunogen Catalog Number: AG25729	Full Name: proline, glutamate and leucine rich protein 1	
	Calculated MW: 1130 aa, 120 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

