

For Research Use Only

IL-4RA/CD124 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67051-2-PBS



Basic Information

Catalog Number:

67051-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28411

GenBank Accession Number:

NM_000418

GeneID (NCBI):

3566

UNIPROT ID:

P24394

Full Name:

interleukin 4 receptor

Calculated MW:

90 kDa

Purification Method:

Protein G purification

CloneNo.:

1E9B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50242-1, IL4R/IL-4RA/CD124 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-3-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50242-4, IL4R/IL-4RA/CD124 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-4-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 0.391-25 ng/mL