

For Research Use Only

Tissue Factor Monoclonal antibody

Catalog Number: 67056-1-Ig **2 Publications**



Basic Information

Catalog Number:

67056-1-Ig

Size:

2422 μ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27563

GenBank Accession Number:

BC011029

GeneID (NCBI):

2152

UNIPROT ID:

P13726

Full Name:

coagulation factor III (thromboplastin, tissue factor)

Calculated MW:

295 aa, 33 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

2E4A1

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : A431 cells, human blood tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Meng-Yun Xu	38886550	Acta Pharmacol Sin	WB
Yiyao Qi	38826705	Heliyon	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 67056-1-Ig (Tissue factor antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.