

For Research Use Only

Tissue Factor/CD142 Monoclonal antibody, PBS Only

Catalog Number: 67056-1-PBS



Basic Information

Catalog Number:

67056-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27563

GenBank Accession Number:

BC011029

GeneID (NCBI):

2152

UNIPROT ID:

P13726

Full Name:

coagulation factor III (thromboplastin, tissue factor)

Calculated MW:

295 aa, 33 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

2E4A1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 67056-1-Ig (Tissue Factor antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67056-1-PBS in a different storage buffer formulation.