

For Research Use Only

REDD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67059-7-PBS



Basic Information

Catalog Number: 67059-7-PBS	GenBank Accession Number: BC007714	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 54541	CloneNo.: 4A10B2
Source: Mouse	UNIPROT ID: Q9NX09	
Isotype: IgG2b	Full Name: DNA-damage-inducible transcript 4	
Immunogen Catalog Number: AG0965	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50649-4, REDD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67059-2-PBS. Detection antibody: 67059-7-PBS. Standard:Ag0965. Range: 0.098-50 ng/mL