

For Research Use Only

C1QA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67063-3-PBS



Basic Information

Catalog Number:

67063-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28622

GenBank Accession Number:

BC030153

GeneID (NCBI):

712

UNIPROT ID:

P02745

Full Name:

complement component 1, q subcomponent, A chain

Calculated MW:

245 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B6D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51151-1, C1QA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67063-2-PBS. Detection antibody: 67063-3-PBS. Standard:Ag28622. Range: 1.563-25 ng/mL