

For Research Use Only

ADRBK1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67067-2-PBS



Basic Information

Catalog Number:

67067-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5071

GenBank Accession Number:

BC037963

GeneID (NCBI):

156

UNIPROT ID:

P25098

Full Name:

adrenergic, beta, receptor kinase 1

Calculated MW:

689 aa, 80 kDa

Purification Method:

Protein G purification

CloneNo.:

2F2A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51612-1, ADRBK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67067-2-PBS. Detection antibody: 67067-1-PBS.



Cytometric bead array standard curve of MP51612-1, ADRBK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67067-2-PBS. Detection antibody: 67067-1-PBS. Standard:Ag5071. Range: 3.125-400 ng/mL