

For Research Use Only

# BCAS3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67085-2-PBS



## Basic Information

Catalog Number: 67085-2-PBS	GenBank Accession Number: BC001250	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 54828	CloneNo.: 1E1D4
Isotype: IgG1	UNIPROT ID: Q9H6U6	
Immunogen Catalog Number: AG28545	Full Name: breast carcinoma amplified sequence 3	
	Calculated MW: 99 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Cytometric bead array sample test of MP51587-1, BCAS3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67085-2-PBS. Detection antibody: 67085-1-PBS.



Cytometric bead array standard curve of MP51587-1, BCAS3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67085-2-PBS. Detection antibody: 67085-1-PBS. Standard:Ag28545. Range: 0.391-50 ng/mL.