

For Research Use Only

RUNX1T1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67086-2-PBS



Basic Information

Catalog Number:

67086-2-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7893

GenBank Accession Number:

BC005850

GeneID (NCBI):

862

UNIPROT ID:

Q06455

Full Name:

runt-related transcription factor 1;
translocated to, 1 (cyclin D-related)

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

1F4C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51598-1, RUNX1T1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67086-2-PBS. Detection antibody: 67086-3-PBS.



Cytometric bead array standard curve of MP51598-1, RUNX1T1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67086-2-PBS. Detection antibody: 67086-3-PBS. Standard: Ag7893. Range: 0.781-100 ng/mL.