

For Research Use Only

# RALB Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67094-2-PBS

## Basic Information

Catalog Number: 67094-2-PBS	GenBank Accession Number: BC018163	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 5899	CloneNo.: 1G7F2
Source: Mouse	UNIPROT ID: P11234	Affinity: $K_D = 1.10 \times 10^{-9}$
Isotype: IgG1	Full Name: v-ral simian leukemia viral oncogene homolog B (ras related; GTP binding protein)	$K_{Off} = 1.89 \times 10^{-5}$ $K_{On} = 1.72 \times 10^4$
Immunogen Catalog Number: AG28600	Calculated MW: 206 aa, 23 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50129-1, RALB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67094-2-PBS. Detection antibody: 67094-1-PBS. Standard:Ag28600. Range: 6.25-100 ng/mL