

For Research Use Only

# ATG9A Monoclonal antibody, PBS Only



Catalog Number: 67096-1-PBS

## Basic Information

Catalog Number: 67096-1-PBS	GenBank Accession Number: BC065534	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 79065	CloneNo.: 1D12E11
Source: Mouse	UNIPROT ID: Q7Z3C6	
Isotype: IgG1	Full Name: ATG9 autophagy related 9 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG24278	Calculated MW: 837 aa, 94 kDa	
	Observed MW: 94 kDa	

## Applications

Tested Applications:  
WB, Indirect ELISA, IHC, IF

Species Specificity:  
Human

## Background Information

ATG9A is the only transmembrane ATG protein essential for autophagy. It plays a key role in the organization of the preautophagosomal structure/phagophore assembly site (PAS). It has been reported that ATG9A expression is increased in oral squamous cell carcinoma and breast cancers. The inhibition of ATG9A can lead to an inhibition of cancer cell proliferation and invasion. (PMID: 29437695, 29568063)

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67096-1-Ig (ATG9A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67096-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using 67096-1-Ig (ATG9A antibody) at dilution of 1:100 and CoraLite488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67096-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using 67096-1-Ig (ATG9A antibody) at dilution of 1:600 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67096-1-PBS in a different storage buffer formulation.