

For Research Use Only

# ELAVL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67097-2-PBS



## Basic Information

Catalog Number:

67097-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28310

GenBank Accession Number:

BC030692

GeneID (NCBI):

1993

UNIPROT ID:

Q12926

Full Name:

ELAV (embryonic lethal, abnormal vision, Drosophila)-like 2 (Hu antigen B)

Calculated MW:

39 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H1A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50972-1, ELAVL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67097-1-PBS. Detection antibody: 67097-2-PBS. Standard:Ag28310. Range: 0.391-100 ng/mL.