

For Research Use Only

PD-ECGF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67104-4-PBS



Basic Information

Catalog Number: 67104-4-PBS	GenBank Accession Number: BC018160	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 1890	CloneNo.: 1F6D1
Source: Mouse	UNIPROT ID: P19971	
Isotype: IgG1	Full Name: thymidine phosphorylase	
Immunogen Catalog Number: AG28631	Calculated MW: 482 aa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51270-2, PD-ECGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67104-4-PBS. Detection antibody: 67104-5-PBS. Standard:Ag28631. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP51270-3, PD-ECGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67104-4-PBS. Detection antibody: 67104-6-PBS. Standard:Ag28631. Range: 1.563-100 ng/mL