

For Research Use Only

PD-ECGF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67104-5-PBS



Basic Information

Catalog Number:

67104-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28631

GenBank Accession Number:

BC018160

GeneID (NCBI):

1890

UNIPROT ID:

P19971

Full Name:

thymidine phosphorylase

Calculated MW:

482 aa, 50 kDa

Purification Method:

Protein G purification

CloneNo.:

1G10G12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51270-2, PD-ECGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67104-4-PBS. Detection antibody: 67104-5-PBS. Standard:Ag28631. Range: 1.563-100 ng/mL.